
How to Update BactoSense Software

Software 2.x.x

1. Insert the update stick (only use the supplied update stick!).
2. Install dotNetDiskImager (Windows only).
3. Launch dotNetDiskImager.
4. Select the image file (*.img).
5. Select the USB device.
6. Check "Verify data when finished."
7. Click "Write to device" and wait until the process is complete.
8. If Windows prompts you to format the USB device, cancel the pop-up without formatting.
9. Eject the USB stick in Windows before removing it.
10. Power off the instrument.
11. Plug the USB stick into the instrument.
12. Power on the instrument.
13. Press "Update" and follow the instructions.

Software 3.x.x

1. Insert any USB stick into your computer (4 GB to 128 GB).
2. Install dotNetDiskImager (Windows only).
3. Launch dotNetDiskImager.
4. Select the image file (*.img).
5. Select the USB device.
6. Check "Verify data when finished."
7. Click "Write to device" and wait until the process is complete.
8. If Windows prompts you to format the USB device, cancel the pop-up without formatting.
9. Eject the USB stick in Windows before removing it.
10. Plug the USB stick into any USB port on the instrument.
11. Navigate to "System parameters" -> "System Update."
12. Wait approximately 10 seconds for the system to mount and decrypt the software.
13. Press "Refresh USB disks."
14. Select the update and press "Update."